

In vitro ADME & PK

PAMPA

Background Information



'The parallel artificial membrane permeability assay (PAMPA), first introduced by Kansy *et al.*, has been widely used in the pharmaceutical industry as a high throughput permeability assay to predict oral absorption.'

¹Di L ,Kerns EH, Fan K, McConnell OJ, and Carter GT. (2003) *Eur J Med Chem* **38**; 223-232.

- The Parallel Artificial Membrane Permeation Assay (PAMPA) is used as an *in vitro* model of passive transcellular permeation.
- PAMPA avoids the complexities of active transport, allowing test compounds to be ranked based on a simple permeability property alone.
- The ability of this assay to evaluate permeability over a large pH range is valuable for an early understanding how new oral compounds might be absorbed across the entire gastrointestinal tract.

Protocol

Test Article Concentration

10 μ M

Number of Replicates

4

Membrane Composition

Hexadecane in hexane (5 % v/v)

Incubation Time

5 hours

Temperature

Room temperature

Test Article Requirements

100 μ L 10 mM DMSO solution

Integrity Marker

Lucifer Yellow

Analysis Method

LC-MS/MS quantification

Data Delivery

P_{app}
Recovery

PAMPA can quickly provide information about passive permeability that is not complicated by other mechanisms such as paracellular transport, active transport and metabolism.



PAMPA

26 Compounds were screened in Cyprotex's PAMPA in quadruplicate on 3 separate occasions. The assay generates consistent, reproducible data over a range of permeability values.

Cyprotex's PAMPA has been successfully trialed by one of our partners and the data compare well with the customer-generated data.

Figure 1

Cyprotex's PAMPA measures passive diffusion of a test compound through an artificial hexadecane membrane.

The protocol was designed in collaboration with our biotechnology partners and follows the method described by researchers at Novartis².




Figure 2

The graph shows the reproducibility of data generated in Cyprotex's PAMPA over 3 separate assays (error bars represent the standard deviation of quadruplicate incubations).





Table 1

Cyprotex $\log P_{app}$ data show a high level of similarity to the third party $\log P_{app}$ data for the purpose of compound classification into low and high $\log P_{app}$.

Compound name	Mean $\log P_{app}$ (Cyprotex)	Mean $\log P_{app}$ (Third Party)
Acyclovir	<-5.79	-6.86
Digoxin	<-6.60	-4.85
Ceftriaxone	<-5.89	-6.25
Fluvastatin	-7.24	-6.31
Ranitidine	-7.16	-6.88
Chloramphenicol	-7.13	-7.76
Amiloride	-7.04	-6.99
Sulfasalazine	-6.70	-6.85
Guanabenz	-5.32	-5.25
Naproxen	-5.04	-5.50
Antipyrine	-4.90	-5.00
Quinidine	-4.55	-4.80
Lansoprazole	-4.45	-4.53
Verapamil	-3.80	-3.59
Desipramine	-3.75	-3.47
Testosterone	-3.69	-3.35

Figure 3

Comparison of Cyprotex PAMPA data with third party data.



Cyprotex data correlate well with data generated by one of our partners ($R^2 = 0.93$). The third party data were produced using the same method as Cyprotex with the exception that a prolonged incubation period and UV absorbance end-point were used. Cyprotex method uses the more sensitive LC-MS/MS end point and a 5 hour incubation.

References

- Di L *et al.* (2003) *Eur J Med Chem* **38**; 223-232.
- Wohnsland F and Faller BJ (2001) *Med Chem* **44**; 923-930.